
Flow-based



**Intervention
Surgical Revision**

**Intragraft Blood Flow
and/or
Intraoperative Blood Flow**

Flow-based AV Access Intervention and/or Surgical Revision

Routine Transonic Flow-QC® surveillance during hemodialysis trends declining vascular access flows to detect the presence of hemodynamically significant stenoses or to identify high flows that threaten cardiac overload.¹⁻³ When a significant stenosis has been identified, various techniques including an endovascular intervention or a surgical revision are used to restore access flow to an acceptable level for hemodialysis delivery. This chapter presents the use of flow measurement with these two well established treatment modalities.⁴

1. Flow-guided Percutaneous Transluminal Angioplasty (PTA)

When AV fistulas or grafts fall below a critical threshold, percutaneous transluminal angioplasty (PTA) is the front line treatment for repairing the access.^{5,6,7,8} Elective rather than emergent, less invasive than surgery, angioplasty conserves the current access site and preserves future access sites. During angioplasty, on-the-spot intragraft flow measurements with the Transonic® ReoCath™ Flow Catheter and Endovascular Flowmeter (HVT100) guides interventional radiologists as they attempt to restore access patency.

2. Flow-based Access Revision Surgery

Surgery is also used to revise challenging AV accesses and restore flow. A number of surgical protocols are used to adjust access flows that can be either too low or too high. Transonic intraoperative flow measurements inform the surgeon and assist in intraoperative decision making during these procedures.

A. Measuring Intra-graft Blood Flow during Angioplasty

After apparently successful angioplasties, access blood flow fails to normalize in one third of patients.¹³⁻¹⁴ Vesely attributes these poor outcomes to the presence of an arterial inflow stenosis, failure to identify multiple lesions within the access, elastic recoil immediately following the angioplasty and/or poor cardiac output. On-the-spot measurements of intra-graft blood flows during the angioplasty procedure (Fig. 3.1) informs the interventionalist about the entire vascular circuit.¹⁵⁻¹⁶

McCarley et al reported that, when surveillance was combined with expedient angioplasty, the rate of graft thrombosis decreased almost fourfold, overall costs dropped by approximately 50% and catheter use decreased dramatically.¹¹ During angioplasty, intra-graft flow measurements provide quantitative data at three points of the intervention: at its outset; during the intervention, and/or at its conclusion.

At the Intervention's Outset

At the procedure's outset, an intra-graft flow measurement confirms, first of all, that the procedure is indeed necessary and blood flows are lower than 600 ml/min K/DOQI-recommended threshold.

Studies have demonstrated that when a static pressure ratio is used for access surveillance, a high flow access can be erroneously targeted for an intervention.¹⁷ It has also been shown that angiographic findings alone do not correlate with blood flow.^{18,19} Vesely addresses this failure of angiography by advising that the intervention should not be performed unless clinical or hemodynamic abnormalities correlate with angiographic findings.⁵

A flow measurement at the beginning of the intervention also provides a flow baseline from which a target flow goal for the procedure can be set. In essence, each individual then serves as his or her own control as the goal of the intervention is improve flow through the access. Intra-graft blood flow measurements at the onset of angioplasty both confirm the need for the intervention and provide a baseline from which flow should improve.




Fig. 3.1: ReoCath™ antegrade flow catheter measuring intra-graft flow before angioplasty.

Flow-based AV Access Intervention *cont.*

During the Intervention

As the intervention proceeds, intragraft ReoCath™ flow measurements provide objective functional data to guide the interventionalist to achieve the target flow goal (see Section C, “How the ReoCath™ System Works,” page 57).

Since Fistula First has mandated increased use of native fistula in lieu of prosthetic grafts, the interventionalist can no longer assume that the cause for declining vascular access flow is a single venous outflow stenosis common in grafts, and correcting it endovascularly will solve the access flow problem. The 2006 K/DOQI Update states, “Inflow stenosis is more common than previously believed. ...Access inflow stenosis occurs in one third of patients referred to interventional facilities with clinical evidence of venous stenosis or thrombosis.”^{7,18,19} Vesely concurs, “Arterial anastomotic stenoses are the single most important flow determinant of the entire vascular access.”⁵

When a Reocath™ intragraft flow measurement demonstrates that the flow target has been attained, the intervention can be safely concluded. However, if, after angioplasty, flow remains below expected levels, the interventionalist is immediately alerted that another problem within the circuit exists that needs to be addressed.

At the Intervention’s Conclusion

Although improvement in blood flow following angioplasty should be obvious, Beathard warns, “Elasticity of venous stenotic lesions appears to be a much more serious problem than has been previously recognized ... Recoil of an elastic lesion can be easily missed.”⁷ To this end, a followup intragraft flow

measurement several minutes after completion of the intervention will confirm continuing good flow and the absence of elastic recoil (Fig. 3.2). Conversely, it will alert the interventionalist if recoil has occurred. Flows can then be documented for the patient record and compared with vascular access flows obtained during subsequent access surveillance in the hemodialysis unit with the Flow-QC® Monitor.




Fig. 3.2: ReoCath™ antegrade flow catheter measuring intragraft flow after angioplasty.

1. Angioplasty Success

Angioplasty success is defined as a significant improvement in access blood flow.^{18,20}

It has been suggested that post-PTA access flow should ideally return to or be higher than former optimum access flow levels, the literature reports the average flow increase to be about 300 ml/min. Krivitski reports that patients whose blood flow increased more than 300 ml/min during angioplasty also showed significant increases in flow measured during subsequent hemodialysis. In his review of 17 studies, the post angioplasty weighted mean flow increase was 309 ml/min for

AV grafts and 329 ml/min for AV fistulas respectively.²² Murray and colleagues also have reported that blood flow achieved after angioplasty is predictive of subsequent graft patency (Fig. 3.3) and asserts that grafts with post-angioplasty access flows of less than 1L/min are more likely to require repeat intervention, to exhibit thrombosis within the first six months, and have a lower one-year survival compared with those with access blood flow rates greater than 1L/min.²⁰

K/DOQI defines successful angioplasty as one where the residual stenosis is less than 30% of the diameter of the access. Ahya et al dispute this anatomic measure, "Visual assessment of the lesion following angioplasty fails to predict the hemodynamic success of the procedure."¹³ A Murray study in which flow actually decreased after angioplasty in four grafts also supports the conclusion that anatomic assessment by angiography does not correlate with flow. Although residual stenoses were visualized in the grafts, none was greater than 30% of the diameter of the access.²¹ Beathard calls for total elimination of the stenosis and cites a Lilly study of 330 cases in which the median primary patency for a graft was longer (6.9 months) when there was no residual stenosis than if there was any degree of residual stenosis (4.6 months).^{7,8}




Fig. 3.3: Murray B et al: "Access Flow after Angioplasty Predicts Subsequent Arteriovenous Graft Survival"¹⁷

Flow-based AV Access Intervention *cont.*

2. How the ReoCath™ System Works

The ReoCath™ Flow Catheter System consists of single-use antegrade and retrograde catheters, a ReoCath™ extension cable and the Transonic® HVT100 Flowmeter (Fig. 3.4). The 6 French antegrade and retrograde flow catheters each have an external injection port connected to a central lumen (Figs. 3.5, 3.6) through which room temperature saline is released into the access during angioplasty.

Each catheter has two temperature sensors (thermistors). When room temperature saline is injected into the access, a thermistor located close to the proximal end of the catheter records the temperature of the injected saline solution. The second thermistor located close to the distal tip of the catheter records the thermodilution within the access.

The catheter connects to the HVT100 Endovascular Flowmeter via a two meter extension cable. The flowmeter automatically calculates and displays intra-graft



Fig 3.4: HVT100 Endovascular Flowmeter measures intragraft blood flow in the arteriovenous (AV) vascular access to provide quantitative information about access functionality during angioplasty.




Fig. 3.5: Antegrade catheter (6 F, 35 cm length) is inserted in the same direction as blood flow. Saline is released proximal to the catheter tip and then is measured downstream by the dilution thermistor.




Fig. 3.6: Retrograde catheter (6 F, 48 cm length) is inserted against the direction of blood flow. Saline is released at the catheter tip and is then is measured downstream by the dilution thermistor.

3. Intra-graft Measurement Protocol

With the ReoCath™ Flow Catheter system, the interventional radiologist introduces a 10cc bolus of room temperature, isotonic saline into the access while performing the endovascular procedure (see Procedure, page 59). The system calculates intra-graft blood flow from the change in temperature of the injected saline. The Flow Measurement Protocol (Flow Charts 3.1, 3.2, pages 60-61) recommends that these intra-graft flows be measured before and after balloon insertion. By comparing the post-angioplasty value to the pre-angioplasty baseline flow value, the interventional radiologist has immediate feedback on the procedure's success. If intra-graft blood flow has not increased to satisfactory levels, the balloon can be re-inserted until flow is optimized.

Principle of Operation

For flow measurements within the AV vascular access, the HVT100 Endovascular Flowmeter and ReoCath™ Flow Catheter use classical dilution-based equations adapted to the unique hemodynamic conditions that exist within the access. Intra-access blood flow measurements are based upon the following equation:

$$Q = k (T_b - T_i) V/S - 0.5 V/\tau$$

Where:

Q = intra-access blood flow;

k = a coefficient related to the thermal properties of blood, saline = 1.08

T_b = temperature of the blood prior to injection;


T_i = temperature of injected saline;

V = volume of injected saline (10ml);

S = the area under the temperature-time dilution curve resulting from the mixing of blood and injected saline;

τ = width of the dilution curve at 50% height (Figure).

The expression (0.5V/τ) is an average expected increase in blood flow as a result of the saline injection.



Thermal dilution curve generated by the change in temperature between the isotonic saline injected into the AV access and the diluted temperature registered by the catheter thermistor within the access.

Flow-based AV Access Intervention *cont.*


INTRA-GRAFT FLOW MEASUREMENTS PROCEDURE

1. Connect HVT100 Endovascular Flowmeter to grounded power receptacle.
2. Select either an Antegrade or Retrograde sterile ReoCath™ Flow Catheter, open pouch and pass the connector to the non-sterile field.
3. Remove the distal cap, curve retainer, and tag from catheter.
4. Have someone outside the sterile field attach the connector of the catheter to the extension cable. Connect extension cable to the HVT100 Endovascular Flowmeter.
Note: The system will not identify the type of catheter (antegrade or retrograde) until it is inserted into the introducer sheath.
5. Open stopcock on the ReoCath™ Flow Catheter and prime the catheter with isotonic saline. Close the stopcock.
6. Insert the catheter through the 6F or larger introducer sheath until the marker band is visible outside the sheath.
7. Press the start button on the front of the Flowmeter.
8. The meter will display “Wait”. The catheter indicator light will display which catheter type is connected, either antegrade or retrograde.
9. After 15-20 seconds, the Flowmeter will display “Ready”. Open the catheter stopcock and inject about 10 ml of room temperature saline (20-25°C) as a 2-3 second bolus in a smooth, continuous motion. (Small deviations in the injection volume do not affect measurement accuracy.)
10. Close the stopcock.
11. The Flowmeter display will change to “PROC. 21” and countdown to “PROC. 00” and will then change to “CALC. 10” (at this point the Flowmeter has begun to calculate flow) The countdown will continue to “CALC 00”. *Note: at some point during the CALC. process, the display may pause momentarily (this is normal).*
12. The Flowmeter will display blood flow in ml/min. If there was a problem with the injection, the display will instead read “REPEAT”. Please return to Step 7.
14. Repeat Steps 7-11 once or twice to ensure measurement reproducibility.
15. Remove the catheter from the introducer sheath and keep it in the sterile field for possible additional measurements later in the procedure.
16. After use, dispose of the catheter according to standard infection control procedures.
17. The extension cable and meter can be cleaned with alcohol wipes.


Note: Catheters can be used repeatedly during the same procedure.

ReoCath™ Flow Measurement Protocol for PTFE Grafts

[Pre-Intervention]



[Post-Intervention]




Flow Chart 3.1: Protocol for Flow-guided intra-graft blood measurement in PTFE grafts


Flow-based AV Access Intervention cont.

ReoCath™ Flow Measurement Protocol for Fistulas

[Pre-Intervention]



[Post-Intervention]



Flow Chart 3.2: Protocol for Flow-guided intra-graft blood measurement in AV fistulas.

4. Validation Studies

Transonic ReoCath™ Flow Catheter technology has been rigorously validated on the bench, *in vivo* and in clinical studies.^{15,16}

4.1. Bench Test

In the bench test model, catheter flow was measured 397 times while water flow in a graft model was varied from 150 - 1,700 ml/min. Measurements between the ReoCath™ Flow Catheter system and the Transonic® HT109 Volume Tubing Flowmeter demonstrated excellent correlation ($r = .98$). More than 60% of the ReoCath™ measurements were within 5% of true flow as measured by the flowmeter and 95% of the measurements were within 15% of the true flow.

4.2. *In Vivo* Testing

In vivo testing in two adult ewes followed. Eleven intragraft blood flow

measurements, obtained with the ReoCath™ Flow Catheter system, were compared to values measured with a Transonic perivascular volume flow-probe. Again, the two measurement methodologies demonstrated excellent correlation ($r = .99$).




Fig. 3.7: Reproducibility of sequential blood flow measurements made by the ReoCath™ Flow Catheter system during the clinical validation study of 24 patients.




Fig. 3.8: Correlation of post-angioplasty blood flow measurements (HD01 Hemodialysis Monitor) with those obtained with the ReoCath™ Flow Catheter system.

Flow-based AV Access Intervention *cont.*

4.3. Clinical Study

In a prospective clinical study, intragraft measurements in 25 patients with the ReoCath™ Flow Catheter system were compared to pre- and post-angiography measurements with the Transonic HD01 Hemodialysis Monitor.²² All patients had PTFE loop grafts, 20 of which were located in the forearm and five were in the upper arm. Fistulograms identified 40 hemodynamically significant stenoses in 24 patients. In each of those patients, angioplasty was performed. Prior to the balloon insertion the ReoCath™ Flow Catheter was inserted through a vascular sheath so that the side holes in the catheter were within the graft, but had not entered its stenotic segment. For each measurement, 10 ml of sterile room temperature, isotonic saline was injected into the catheter's injection port over 3-4 seconds. Two consecutive intragraft blood flow measurements were taken and the results averaged. If the two measurements differed by more than 10%, a third measurement was taken. Sequential measurements were highly reproducible (Fig. 3.7). Measurements were performed both before and immediately after angioplasty.

On an average of 11.9 days before angioplasty, the mean vascular access blood flow measured by the HD01 Monitor in the 24 subjects was 463 ± 154 ml/min. When intragraft flow was measured just prior to balloon insertion with the ReoCath™ Flow Catheter system, the mean flow was 495 ± 180 ml/min. Post-angioplasty correlation between the ReoCath™ flows and HD01 flow measurements were 779 ± 331 ml/min and 781 ± 221 ml/min respectively (Fig. 3.8). The mean increase in vascular access blood flow following angioplasty measured by the ReoCath™ system was 324 ± 267 ml/min and 319 ± 256 ml/min measured by the HD01 Monitor 4.9 days (average) following angioplasty.

4.4. Conclusion

The ReoCath™ system provides a quantitative flow assessment of the entire vascular circuit: its arterial inflow segment, mid-access segment and venous outflow segment.¹⁵⁻¹⁶ This data alerts the interventional radiologist to lesions wherever they are located within the circuit to ensure success of the intervention and avert an immediate return of the patient to the radiology suite or to surgery for another intervention.²³⁻²⁴

Flow-based AV Access Revision Surgery

B. Flow-based AV Access Revision Surgery

Creation of AV fistulas or grafts produces unique cardiovascular hemodynamics. Oxygenated arterial blood usually flows through an intricate high resistance system of arterioles and capillaries before it reaches the venous system. When an AV access is created and arterial blood flows through the access, it immediately encounters a low resistance vein. This makes AV accesses (prosthetic grafts, in particular) susceptible to myointimal hyperplasia that can eventually lead to stenosis, thrombosis and access failure.

A well functioning access remains free of complications and provides sufficient flow to deliver the prescribed dialysis prescription. When access flow is inadequate, options include: repairing it endovascularly, and/or surgically revising the access, or a creating a new access at another site. When surgery is chosen, intraoperative flow measurements during the revision provide quantitative flow data to help the procedure attain its desired conclusion.

“Accurate flow measurements can be of great assistance during vascular reconstructive surgery. The primary aim with these intraoperative measurements is to obtain information on the immediate result of the reconstruction, where a technical failure may jeopardize an otherwise successful operation.”²⁵

1. Low Access Flow

If an access fails or a patient develops acute upper arm or extremity ischemia, the vascular surgeon must repair the access or establish another access rather than resort to use of a catheter for long-term dialysis. During these surgeries, intraoperative flow measurements help guide the revision to a successful outcome.

When the more common AV access sites are exhausted, vascular surgeons may turn to creative complex bypasses in order to establish a viable AV access for administration of hemodialysis.²⁸⁻³⁰

Examples of complex expanded polytetrafluoroethylene (ePTFE) bypass grafts reported in the literature include: a necklace bypass (axillary artery to contralateral axillary vein); contralateral internal jugular vein (brachial artery to the internal jugular vein); femorofemoral crossover or “Bikini Bypass” (femoral artery to contralateral femoral vein); axillary artery to popliteal vein (axillary artery to superficial femoro-popliteal vein); femoral artery to right atrium.²⁸⁻³⁰

Flow-based AV Access Revision Surgery *cont.*

When the ESRD patient returns to hemodialysis, ongoing quantitative surveillance of a revised access can then be resumed with ultrasound dilution monitoring of graft inflow rate, recirculation, access flow and cardiac output described in Chapter II of the handbook.

The following two case studies demonstrate the use of intraoperative flow measurements as a surgeon revises problematic accesses with creative solutions.

1.1: Case Study

Brachiocephalic Fistula with Cephalic Arch Stenosis: Correction by Cephalic Vein Ligation and Anastomosis to the Axillary Vein

A 58-year-old male patient, on dialysis for the past 25 years, has had three transplants that had ultimately failed. He had a very high rate of panel reactive antibodies (PRA) which made him unsuitable for another transplant.

Attempts at peritoneal dialysis (PD) ended in peritonitis. He was being dialyzed through a precious 12-year-old right brachiocephalic arteriovenous fistula (AVF). The patient was well aware of his disease and very cooperative.

Cephalic Arch Stenosis

During the previous year, his vascular access flow measurements, as measured by the Transonic Flow-QC® Hemodialysis Monitor, dropped. He was referred for a fistulogram that revealed a right cephalic arch stenosis (Fig. 3.9). The stenosis was successfully dilated three times, but, the intervals between the dilations and the next access flow decrease was only eight weeks.

Clinically, the fistula was pulsatile with no thrill or bruit. Dialysis was painful and inefficient with very high venous pressures, low pump speed and poor inflow (200-250 ml/min). It was decided to correct the problem surgically through creation of a new arch, by ligating the cephalic vein, just proximal to the stenosis, mobilizing it and swinging it over to anastomose with the axillary vein.




Fig. 3.9: Fistulogram showing cephalic arch stenosis.




Fig. 3.10: Right brachiocephalic fistula with cephalic arch stenosis.

continued on next page

2. High Access Flow

2.1: Cardiac Overload

The unique hemodynamics produced by AV fistula creation may cause AV fistulas to mature to flow levels (greater than 2L/min) that can overwork the heart and lead to congestive heart failure. When AV access flows exceed 2 liters per minute, cardiac output measurements with the Flow-QC® Hemodialysis Monitor and Central Hemodynamic Profiling inform whether the heart is indeed stressed (see pages 44-48).

Banding

When cardiac overload occurs, the traditional approach has been to band or lengthen the fistula to increase resistance and reduce cardiac overload. When the surgeon chooses to “band” a fistula, intraoperative flow measurements confirm that “banding” has reduced flow to an acceptable level.³¹⁻³²

Inflow Reduction

Banding of a fistula, however, can be counterproductive. Although fistula flow is reduced and less stress is placed on the heart, the reduction in blood flow within the access might accelerate stenosis and subsequent thrombosis. Recognition of this danger has led to the use of alternative methods to reduce flow through an high access bypass.

A 2007 study by Chemla and colleagues reports anastomotic distalization to reduce inflow in the vascular access to treat high cardiac output in hemodialysis patients.³³

The access inflow reduction procedure described in the study, calls for the fistula anastomosis to be exposed, dissected and tied off. The brachial artery is then reconstructed with either an end-to end anastomosis or a patch. The radial artery is then dissected at the wrist and a 6 mm expanded polytetrafluoroethylene (ePTFE) graft is implanted between the radial artery and the cephalic vein and bypass outflow is then measured (Figs. 3.14, 3.15).




Fig. 3.14: Fistula flow is measured.

Flow-based AV Access Revision Surgery *cont.*

The anastomotic distalization decreased access mean flow significantly with a resultant decrease in cardiac output and resolution of symptoms (Table 3.3).

INFLOW REDUCTION SURGERY (n = 17)	MEAN FLOW ML/MIN	CARDIAC OUTPUT L/MIN
PRE-SURGERY	3135 ± 692	8 ± 3.1
POST-SURGERY	1025 ± 551	5.6 ± 1.7

Table 3.3: Mean blood flow and cardiac output in 17 patients pre- and post-inflow reduction surgery.




Fig. 3.15: Following construction of a bypass graft from the radial artery to the cephalic vein, outflow is measured.

“Flow reduction using intraoperative access flow monitoring is an effective and durable technique allowing for the correction of distal ischemia and cardiac insufficiency in patients with a high-flow autogenous access.”³⁴

2.2: Measurement Steps during Intraoperative Fistula Banding or Inflow Reduction Surgery

1) Establish Flow Goal

Before surgery, an AV access flow goal is established, with consideration to the type and location of the access, patient size, weight, gender, medical and access flow histories.

2) Measure Flow before Revision

Measure flow through the AV access before the AV access revision per protocols on pages 11 and 12.

3) Measure Flow after Revision and Compare Flows

Measure access outflows after revision. Compare and document the before and after flows. Evaluate the flow measurement after repair or revision against the desired AV access flow target.

3. Ischemic Steal Syndrome

Ischemic steal syndrome (ISS) following hemodialysis is a devastating complication that occurs when an AV access that shunts arterial inflow into the low-pressure venous circulation “steals” arterial flow from lower arm and hand perfusion. While some level of “steal” is an almost universal “physiologic” result of vascular access creation, it is usually asymptomatic. However, ischemic steal syndrome develops in 1.6 - 8% of patients whose inherent compensatory mechanisms such as collateral circulation and vasodilatation cannot meet metabolic demands.

Ischemic steal syndrome challenges the clinician to relieve distal ischemia while maintaining a functional hemodialysis access. Traditional approaches include banding or lengthening the fistula to increase fistula resistance, but the resulting decrease in blood flow may also lead to access stenosis or thrombosis.

3.1: Alternative Techniques to Relieve Ischemic Steal Syndrome

Distal revascularization-interval ligation (DRIL) technique has become an accepted alternative technique for managing ISS. The procedure eliminates the potential pathway for steal syndrome by ligating the artery distal to the AV fistula origin, and revascularizes the extremity through creation of a bypass (saphenous vein, PTFE graft) from above the AV fistula to below the AV fistula (Fig. 3.16).³⁵⁻⁵⁰ However, because the DRIL technique calls for ligating the artery, it has not become universally accepted.

Wayne Gradman and colleagues from Los Angeles' Cedars Sinai Medical Center used mathematical models and intraoperative flow measurements to investigate several options for mitigating ischemic steal syndrome. He and other have proposed alternative procedures for alleviating ischemic steal syndrome.⁴⁸⁻⁵⁰ These include:




Fig. 3.16: DRIL includes distal artery ligation and bypass construction to an extremity.

Flow-based AV Access Revision Surgery *cont.*

Revision Using Distal Inflow (RUDI)⁵⁰

The fistula is ligated at its origin and then re-established via a bypass from a more distal arterial source to the venous limb of the fistula.

Proximalization of Arterial Inflow (PAI)⁵⁰⁻⁵³

The fistula is ligated at its origin and then re-established via a bypass from a more distal arterial source to the venous limb of the fistula.

3.2: Flow Measurement Steps for Procedures to Relieve ISS

1) Expose Arterial Segment for Pre-Bypass Flow

Expose arterial inflow segment distal to the fistula origin. Select an optimum site (wide enough to accommodate the probe's acoustic reflector) for applying the probe, and clean the site of fat and excess tissue.

2) Select Flowprobe Sizes

Estimate the diameter of vessel with a gauge. Select a probe size so that the vessel diameter will fill 60% - 100% of the flowprobe window.

<u>SITE</u>	<u>PROBE SIZE</u>	<u>NONRESTRICTIVE VESSEL RANGE</u>
RADIAL	2 mm	1.5 - 2.7 mm
ARTERY	3 mm	2.5 - 7.7 mm

3) Apply Flowprobe and Measure Flow

Apply sterile Aquasonic Gel 100™ to the Flowprobe to provide ultrasound coupling between the probe body and probe reflector. Apply the Flowprobe to the site bending the probe's flexible neck segment, as necessary, so that the entire vessel lies within the flowprobe window and aligns with the probe body. Listen to the pitch of FlowSound™ as the Flowprobe is applied to the vessel. The higher the pitch, the greater the flow. Note if flow is going antegrade or retrograde. Retrograde flow indicates ISS.

Check the Signal Quality Indicator (bucket display) on the flowmeter's front panel for ultrasound acoustic contact. An acoustic error message will be displayed if ultrasound contact falls below an acceptable minimum.

4) Construct the Bypass

5) Measure and Evaluate Bypass Flow

With the Flowprobe positioned as under Step 3, measure bypass flow.

3.3: Case Study

ISS on a High Inflow Brachiocephalic Fistula: Management by Inflow Reduction Using a Bovine Bypass Graft from the Radial Artery⁵⁴⁻⁵⁵

Introduction

A female patient with a mature brachiocephalic fistula presented with ischemic steal syndrome (ISS) as she awaited a second nephrectomy. The fistula had not been cannulated. Surgery was undertaken to alleviate ISS symptoms before the nephrectomy.

Baseline Flow Measurement

Fistula (cephalic vein) flow was measured with a Transonic perivascular flowprobe (Fig. 3.17). Flow measured 1600 ml/min - 1700 ml/min confirming a mature fistula.

Fistula Anastomosis Ligation

The anastomosis was ligated and the brachial artery repaired (Fig. 3.18a).

Bypass Graft Construction

A bovine graft was constructed between the radial artery and the cephalic vein (Fig. 3.18b).

Post-Bypass Flow Measurements

Cephalic vein flow after revascularization dropped to 600 ml/min (Fig. 3.18c). Bypass flow also measured approximately 600 ml/min (Fig. 3.18d).

Conclusion: After surgery, steal symptoms immediately began to improve as the brachial artery resumed full delivery of flow to the distal arm and hand circulation (Table. 3.4).




Fig. 3.17: Fistula flow is measured pre revision.




Fig. 3.18: After construction of a radial artery-cephalic vein graft, flow was remeasured.

FLOW SUMMARY		
CONDUIT	BASELINE	POST REVASCLARIZATION
CEPHALIC VEIN	1600 -1700 ml/min	600 ml/min
BYPASS GRAFT		~600 ml/min

Table 3.4: Cephalic vein flow drop from 1600-1700 ml/min to 600 ml/min following inflow reduction surgery.

Case study courtesy of **Eric S. Chemla, MD**, consultant surgeon and honorary senior lecturer, Renal Transplant and Vascular Surgery, St. George's Healthcare NHS Trust, London, UK

C. References

- 1 Beathard, GA, "The Treatment of a Vascular Access Graft Dysfunction: A Nephrologist's View and Experience," *Adv Ren Replace Ther* 1994; 1: 131-147.
- 2 Schwab, SJ et al, "Hemodialysis Arteriovenous Access: Detection of Stenosis and Response to Treatment by Vascular Access Blood Flow," *Kid Int* 2001; 59(1): 358-362. (HD172A)
- 3 Besarab, A, "Access Monitoring Is Worthwhile and Valuable," *Blood Purification* 2006; 24: 77-89.
- 4 Tessitore, N et al, "Endovascular versus Surgical Preemptive Repair of Forearm Arteriovenous Fistula Juxta-Anastomotic Stenosis: Analysis of Data Collected Prospectively from 1999 to 2004," *Clin J Am Soc Nephrol* 2006; 1: 448-454. (HD405 A)
- 5 Vesely, TM, "Is Percutaneous Transluminal Angioplasty an Effective Intervention for Arteriovenous Graft Stenosis, Opinion?" *Seminars in Dialysis* 2005; 18(3): 197-198.
- 6 Asif, A et al, "Percutaneous Management of Perianastomotic Stenosis in Arteriovenous Fistulae: Results of a Prospective Study," *Kidney Int* 2006; 69(10): 1904-9.
- 7 Beathard, GA, "Is Percutaneous Transluminal Angioplasty an Effective Intervention for Arteriovenous Graft Stenosis, Opinion?" *Seminars in Dialysis* 2005; 18(3): 194-196.
- 8 Lilly, RZ et al, "Predictors of Arteriovenous Graft Patency after Radiologic Intervention in Hemodialysis Patients," *Am J Kid Dis* 2001; 37:945-954.
9. Vesely, TM, "Endovascular Procedures: New Techniques and New Technology," *Vascular Access for Hemodialysis—IX*, edited by Mitchell L. Henry, MD, WL Gore & Associates, Inc., Bonus Books, 2005; 13-23, Los Angeles, CA. (HD7067A)
- 10 Tessitore, N et al, "Can Blood Flow Surveillance and Pre-emptive Repair of Subclinical Stenosis Prolong the Useful Life of Arteriovenous Fistulae? A Randomized Controlled Study," *Nephrol Dial Transplant* 2004; 19: 2325-2333. (HD405 A)
- 11 McCarley, P et al, "Vascular Access Blood Flow Monitoring Reduces Access Morbidity and Costs," *Kidney Int* 2001; 60(3): 1164-1172. (HD317A)
- 12 Sands, MS, "Utility of Angioflow Catheters," *VEITH Symposium* 2006.
- 13 Ahya, SN et al, "Flow in Hemodialysis Grafts After Angioplasty: Do Radiologic Criteria Predict Success?" *Kidney Int* 2001; 59(5): 1974-1978. (HD175A)
- 14 Van Der Linden, J et al, "Short- and Long-Term Functional Effects of Percutaneous Transluminal Angioplasty in Hemodialysis Vascular Access," *J Am Soc Nephrol* 2002; 13(3): 715-720. (HD256A)
- 15 Vesely, TM et al, "Use of a Catheter-based System to Measure Blood Flow in Hemodialysis Grafts during Angioplasty Procedures," *J Vasc Interv Radiol* 2002; 13(4): 371-378. (IR1V)
- 16 Vesely, TM et al, "Preliminary Experiences Using Intravascular Blood Flow Monitor (IBFM) during Vascular Access Angioplasty," *J Am Soc of Nephrol Abstracts* 1999; 10: 221A.
- 17 Spergel, LM et al, "Static Intra-access Pressure Ratio Does Not Correlate with Access Blood Flow," *Kidney Int* 2004; 66(4): 1512-1516.
- 18 National Kidney Foundation. *K/DOQI Clinical Practice Guidelines for Vascular Access*. 2006 Update <http://www.kidney.or-/professionals/kdoqi/guideline-upHD-PD-VA/va-guide4.htm>.

Flow-based AV Access Revision Surgery *cont.*

- 19 Asif, A et al, "Inflow Stenosis in Arteriovenous Fistulas and Grafts: A Multicenter, Prospective Study," *Kidney Int* 2005; 67:1986-1992.
- 20 Murray, BM et al, "Access Flow after Angioplasty Predicts Subsequent Arteriovenous Graft Survival," *J Vasc Interv Radiol* 2006; 17: 303-309. (HD2020A)
- 21 Tessitore, N et al, "A Prospective Controlled Trial on the Effect of Percutaneous Transluminal Angioplasty on Functioning Arteriovenous Fistulae Survival," *J Am Soc Nephrol* 2003; 14(6): 1623-1627. (HD310A)
- 22 Krivitski, NM "Access Flow Measurement during Surveillance and Percutaneous Transluminal Angioplasty Intervention," *Seminars Dialysis* 2003; 16(4): 304-308.
- 23 Levine, MI et al, "The AngioFlow Catheter (TC) (Transonic, Inc.) Predicts Increases in AV Access (AVA) Blood Flow (Qa) Post (p) Angioplasty (A) Measured Using the Transonic HD01 Flow Monitor (TM)," *J Am Soc Nephrol Abstracts* 2001; 12: 294A. (HD223A)
- 24 Levine, MI et al, "Patency Rates and Intra-Access Blood Flow (Qa) When Using the Angioflow Catheter (AC) to Verify Arteriovenous Dialysis Vascular Access Percutaneous Angioplasty Success," *J Am Soc Nephrol Abstracts* 2003; 14: 724A. (HD319A)
- 25 Lundell, A, Begqvist, D, "Intraoperative Flow Measurements in Vascular Reconstruction," *Annales Chirurgiae Gynaecologiae* 1992; 81(2): 187-191. (3G)
- 26 Calder, FR et al, "The Axillary Artery-Popliteal Vein Extended Polytetrafluoro-ethylene Graft: A New Technique for The Complicated Dialysis Access Patient," *Nephrol Dial Transplant* 2004; 19(4) 998-1000.
- 27 Chemla, E, et al, "Complex Bypasses and Fistulas for Difficult Hemodialysis Access: A Prospective Single-Center Experience," *Seminars in Dialysis* 2006; 19(3): 246-250.
- 28 Chemla, E, Exotic Surgery 5th International Congress of the Vascular Access Society (VAS). 2007, June 11-13, Nice, France. Abstract L-040 (7471AH)
- 31 Murray, BM et al, "Effect of Surgical Banding of a High-Flow Fistula on Access Flow and Cardiac Output: Intraoperative and Long-Term Measurements" *Am J Kid Dis*, 2004; 44(6): 1090-1096. (HD415A)
- 32 Thermann F et al, "Dialysis Shunt-associated Steal Syndrome (DASS) Following Brachial Accesses: The Value of Fistula Banding under Blood Flow Control," *Langenbecks Arch Surg* 2007. (7473R)
- 33 Chemla, E, et al, "Inflow Reduction by Distalization of Anastomosis Treats Efficiently High-Inflow Cardiac Output Vascular Access for Hemodialysis," *Sem Dialysis* 2007; 20(1): 68-72.
- 34 Zanow, J et al, "Flow Reduction in the High-flow Arteriovenous Access Using Intraoperative Flow Monitoring," *J Vasc Surg* 2006;44: 1273-1278. (7410AHM)
- 35 Mills, JL et al, "Management of Hand Ischemia Associated with Arteriovenous Hemodialysis Access," Chapter 88, *Mastery of Vascular and Endovascular Surgery*. Zelenock, GB, Huber, TS, Messina, LM, Lumsden, AB, Moneta, GL, editors, Lippincott Williams & Wilkins, Philadelphia, 2006.
- 36 Schanzer, H et al, "Treatment of Ischemia Due to "Steal by Arteriovenous Fistula with Distal Artery Ligation and Revascularization," *J Vasc Surg* 1988; 7(6): 770-3.
- 37 Schanzer, H et al, "Treatment of Angio-access-induced Ischemia by Revascularization," *J Vasc Surg* 1992; 16(6): 861-4.
- 38 Berman, SS et al, "Distal Revascularization-Interval Ligation for Limb Salvage and Maintenance of Dialysis Access in Ischemic Steal Syndrome." *J Vasc Surg* 1997; 26(3):93-402; discussion 402-4.

Flow-based AV Access Revision Surgery *cont.*

- 39 Wixon, CL et al, "Distal Revascularization-Interval Ligation for Maintenance of Dialysis Access and Restoration of Distal Perfusion in Ischemic Steal Syndrome," *Semin Vasc Surg* 2000 ;13(1): 77-82.
- 40 Jean-Baptiste, RS, Gahtan V, "Distal Revascularization-interval Ligation (DRIL) Procedure for Ischemic Steal Syndrome (ISS) after Arteriovenous Fistula Placement," *Surg Technol Int* 2004;12: 201-5.
- 41 Knox, RC et al, "Distal Revascularization-Interval Ligation: a Durable and Effective Treatment for Ischemic Steal Syndrome after Hemodialysis Access," *J Vasc Surg* 2002; 36(2):250-5.
- 42 Ascitutto, G et al, "Distal Revascularization-Interval Ligation for the Treatment of Angioaccess-induced Ischemia. Case Report," *Minerva Urol Nefrol*, 2006 Mar;58(1):91-5.
- 43 Diehl, L et al, "Operative Management of Distal Ischemia Complicating Upper Extremity Dialysis Access," *Am J Surg* 2003;186(1): 17-9.
- 44 Sessa, C et al, "Treatment of Hand Ischemia Following Angioaccess Surgery Using the DRIL Technique with Preservation of Vascular Access: Description of An 18 Case Series," *Ann Vasc Surg* 2004;18(6): 685-94.
- 45 Walz, P et al, "Distal Revascularization and Interval Ligation (DRIL) Procedure for the Treatment of Ischemic Steal Syndrome after Arm Arteriovenous Fistula." *Ann Vasc Surg* 2007.
- 46 Illig, KA et al., "Ischemia: DRIL," 5th International Congress of the Vasc Access Soc 2007 Nice France, Extended Abstract L-075 (7466AHR)
- 47 Illig, KA et al., "Hemodynamics of Distal Revascularization Interval Ligation," *Ann Vasc Surg* 2005; 19:199-207. (7228AH)
- 48 Gradman, WS, Pozrikidis C, "Analysis of Options for Mitigating Hemodialysis Access-Related Ischemic Steal Phenomena," *Ann Vasc Surg* 2004; 18: 59-65. (2938AH)
- 49 Tynan-Cuisinier, GS et al, "Strategies for Predicting and Treating Access Induced Ischemic Steal Syndrome," *Eur J Vasc Endovasc Surg* 2006; 32(3): 309-15.
- 50 Minion, DJ et al, "Revision Using Distal Inflow: A Novel Approach to Dialysis-associated Steal Syndrome." *Ann Vasc Surg* 2005; 19(5): 625-8.
- 51 Zanow, J et al, "Proximalization of the Arterial Inflow: A New Technique to Treat Access-related Ischemia," *J Vasc Surg* 2006; 44(5): 1134.
- 52 Zanow, J et al, "Proximalization of the Arterial Inflow to Treat Access-related Ischemia" 5th International Congress of the Vascular Access Society (VAS) 2007, Nice, France, Abstract L-076 (7467AHR)
- 53 Themann, F, Wollert, U, "Proximalization of the Arterial Inflow (PAI)" : New Therapeutic Option in Case of Advanced Steal Syndrome after Brachial Dialysis Access," 5th International Congress of the Vascular Access Society (VAS). 2007, Nice, France, Abstract L-077 (7468AHR)
- 54 Chemla, E et al, "Intraoperative Flow Measurements Are Helpful in the Treatment of High-Inflow Steal Syndrome on a Predialysis Patient with a Brachiocephalic Fistula: A Case Report," *Ann Vasc Surg* 2007.
- 55 Chemla, E, "Interim Results on a Prospective, Randomized Comparison of PTFE Graft Versus Decellularised Bovine Ureter for AV Haemodialysis Access," 5th International Congress of the Vascular Access Society 2007, Nice, France, Abstract L-076 (7465R)

Acknowledgement

Transonic Systems' technological flow innovations include: intraoperative flowmetry using transit-time ultrasound, vascular access flow measurement during hemodialysis via blood line reversal, and catheter-based intragraft flowmetry during AV access intervention. Each of these technologies has spawned significant improvements in hemodialysis patient care.

The three measurement techniques were developed under National Institutes of Health (NIH) grants. We gratefully acknowledge this paramount financial assistance.

The true innovations, however, come from end users: innovation-minded clinicians who are dedicated to the study and exploration of flow-based innovations that improve patient care. We, and the larger medical community, owe much to these innovators. Without their continuing work, the ultimate goal of all ESRD studies – to remove “end-stage” from ESRD – would remain a pipe dream.

We especially thank and acknowledge Eric S. Chemla, MD, consultant surgeon and honorary senior lecturer, Renal Transplant and Vascular Surgery, St. George's Healthcare NHS Trust, London, UK, who took the time to read the manuscript, offer corrections and share cases from his vascular surgical practice.



USA: Tel: 800-353-3569; Fax 607-257-7256; www.transonic.com
Europe: Tel: 31 43 407 7200; email: info@transonic.nl; Fax: 3143 407 7201
Asia: Tel: 886 3399-5806; Fax: 886 3399-5805