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Transonic Hemodialysis Flow-QC® Set: Instructions for Use

Precautions

- The Transonic Flow-QC® Set for Hemodialysis is for single use only.
- Discard Flow-QC® Set as biohazard waste per local or institutional guidelines after removal from patient or after 12 hours of use.
- Reusing a single-use medical device, or using it beyond the prescribed time may be unsafe. Risks include, but are not limited to, cross-infection, contamination, mechanical failure, patient injury or other associated patient health risks.
- Transonic Systems Inc. does not assume any legal liability for the reuse of a single-use medical device or the use of the device beyond the prescribed time.
- Use standard institutional aseptic techniques.
- The device contains Pthalates which are below acceptable limits. The device is acceptable to be used on all patient populations.

Setup Procedures


- Maintaining aseptic technique throughout, place the Flow-QC® Set in the hemodialysis circuit during priming of the dialysis lines with normal saline as follows:
- Open sterile Flow-QC® Tubing set
 - For ARTERIAL connection: Remove the Recirculation Cap from the male luer of the arterial HD blood line. Then remove the red cap from female luer of Flow-QC® tubing set (Fig 1, B) and connect to the arterial HD tubing line.
 - For VENOUS connection: Remove the Recirculation Cap from the male luer of the Venous HD blood line. Then remove the blue cap from the female luer of Flow-QC® tubing (Fig 1, C) and connect to the Venous HD blood line.
- Remove clear plastic caps from both Arterial and Venous sides of Flow-QC® Tubing set (Fig 1, A) and connect the Arterial side to the priming bag and the Venous side to the waste bag.
- Prime entire HD Blood line as usual including Flow-QC® set. Take care to remove all air bubbles from the Flow-QC® Set during priming.
- To ensure that the injection port (Fig 1, D) is primed, open clamp and loosen cap (Fig 1, E) until all bubbles are cleared. When clear, re-clamp line.
- Attach the arterial and venous Flow-QC® Set lines to the fistula needles of the patient as usual at point A.

Delivered Flow, Access Recirculation and Access Flow Measurements:

- Place the venous and arterial flow/dilution sensors on their respective Flow-QC® Set lines (Figure 2 #F).
- See Quick Reference Guide. NOTE: Saline bolus for the Recirculation or Access Flow measurements must be administered via the saline bag release or the venous bubble trap NOT via the Flow-QC® injection port
- Reverse bloodlines only at ends #A going to Fistula Needles (Fig 2).

Cardiac Output Measurements

- Cardiac Output Measurements can only be performed if there is no Access Recirculation.
- Unscrew the cap (#E) on the injection port and attach a needleless syringe filled with 30 mL of normal saline warmed to body temperature (37°C) to the injection port (Fig 3, #D) via the luer-loc connection. DO NOT USE NEEDLES.
- See Quick Reference Guides.
- Make measurement injection quickly and smoothly (~6-8 sec).
- When all measurement injections are complete, place cap (#E) on injection port.
- Flow-QC® tubing can remain in the line for the duration of treatment; discard AFTER removal of dialysis lines from patient.

transonic
THE MEASURE OF BETTER RESULTS.

⚠ READ INSTRUCTIONS IN OPERATOR'S MANUAL

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