## Hemodialysis

## Introducing: **Isopropyl Alcohol Lubricant for** Transonic® Flow/dilution Sensors

- Transonic has validated and approved the use of 70% isopropyl alcohol as a lubricant for Transonic Hemodialysis Flow/dilution Sensors
- Meets system specifications for signal strength and accuracy versus Vaseline®
- Faster, easier procedure with the antiseptic benefits of 70% isopropyl alcohol







Transonic Systems Inc. Tel: +1 607-257-5300 Fax: +1 607-257-7256 support@transonic.com





Transonic Europe B.V. Tel: +31 43-407-7200 Fax: +31 43-407-7201 europe@transonic.com



## Instructions For Use

## 70% Isopropyl Alcohol Lubricant for Transonic Hemodialysis Flowsensors

- 1. Open the door of the first paired Flow/dilution Sensor.
- 2. Place the tubing segment to be inserted next to the Flow/dilution Sensor. The arrow on the Sensor must point in the direction of flow.
- 3. Open a 70% isopropyl alcohol wipe (prep pad).
- 4. Wipe the entire circumference of the tubing segment which will be inserted into the Flow/dilution Sensor.
- 5. Immediately insert this tubing segment into the Flow/dilution Sensor and close the door. The alcohol evaporates from the tubing surface in a few seconds. It is important to wipe, insert tubing, and close door as a continuous sequence.
- 6. Repeat for the second paired Flow/dilution Sensor and tubing segment.
- 7. When the Hemodialysis Monitor has been turned on and a patient is being measured, a Signal Strength indicator will display in the upper left of the Monitor screen. Verify that this Signal Strength indicator is green. This means that the paired Flow/dilution Sensors have adequate contact with the tubing. If the Signal Strength indicator is not green, repeat the [wipe-insert-close] sequence to achieve proper contact.