

# Optimizing Dialysis Adequacy In Dialysis Patients with Central Venous Catheters


Even though central venous catheters (CVCs) are prone to thrombosis and infection, 80.9% of patients use a catheter at initiation of HD, and 21.1% of prevalent patients continue that use.<sup>1</sup> (Note the 2019 KDOQI updated definition of CVC dysfunction: failure to maintain the prescribed extracorporeal blood flow required for adequate hemodialysis without lengthening the prescribed HD treatment.) Yet, catheter dysfunction remains a serious cause for concern for hemodialysis providers. Two potential pitfalls to achieving adequate catheter dose delivery include:

- A fibrin sheath can block the catheter's lumen, thus impeding flow and causing a severe drop in dialysis dose delivery.
- The close proximity of the catheter's arterial entry and venous return ports make recirculation and underdialysis likely.

**Note:** if any intervention occurs, such as the use of a thrombotic agent, the Delivered Flow and Recirculation measurements can be repeated to determine the effectiveness of the intervention.

## Delivered Blood Flow and Recirculation to Optimize Catheter Dialysis Measurements

Compare Transonic Delivered Blood Flow reading with the hemodialysis machine's pump setting. If the disparity is more than 10%, check for kinked tubing. A fibrin sheath might be restricting inflow and reducing dose delivery. The optimization of hemodialysis for catheter connection configuration with the Transonic Hemodialysis Monitor can be used to then check for recirculation. If the connection is then reversed, the Delivered Flow and Recirculation measurements should then



Central venous catheter inserted via the jugular vein into the right atrium of the heart to serve as a vascular access for hemodialysis.

be repeated to determine the best catheter configuration.

- The nurse can adjust the dialysis delivery parameters (time, pump setting etc.) to compensate for recirculation and deliver the prescribed dose of dialysis to the patient.
- Dialysis lines may be reversed. Reversing the lines might also correct high recirculation.

The nurse should report unusual delivered blood flow and recirculation readings to the Patient Care Team and/or nephrologist to ensure optimum short- and long-term management of the patient's hemodialysis treatment.

1. United State Renal Data System. 2022 USRDS Annual Data Report: Epidemiology of kidney disease in the United States. National Institute of Health, National Institute of Diabetes and Digestive and Kidney Diseases. Bethesda, MD, 2022.

# Optimizing HD Adequacy in Catheters

## Catheter Configuration with the Transonic HD Monitor

### Step 1:

#### MEASURE DELIVERED BLOOD FLOW RATE

With the bloodlines configured as normally used (document configuration), measure flow. Transonic Delivered Blood Flow Rate (Qb) is within 0-10% of the hemodialysis machine's set blood pump speed or delivery flow rate.\*

YES

TRANSONIC DELIVERED BLOOD FLOW RATE (QB) IS WITHIN 0-10% OF HEMODIALYSIS MACHINE'S SET, OR DELIVERY FLOW READING\*

Current blood pump setting is maximizing the Delivered Blood Flow with the current catheter to bloodline configuration.

PROCEED TO RECIRCULATION MEASUREMENT

YES

TRANSONIC DELIVERED BLOOD FLOW RATE (QB) IS WITHIN 0-10% OF THE HEMODIALYSIS MACHINE'S SET OR DELIVERY FLOW READING\*

Current blood pump setting is maximizing the Delivered Blood Flow with the current catheter to bloodline configuration.

PROCEED TO RECIRCULATION MEASUREMENT

\*Some Hemodialysis machines display both a Set Blood Pump Speed and Delivery Flow Reading. If both readings are displayed on your Hemodialysis machine, use the Delivery Flow Reading.

#### Catheter Configurations:

- Normal Configuration: Arterial Catheter Hub to Arterial Bloodline + Venous Catheter Hub to Venous Bloodline
- Reverse Configuration: Arterial Catheter Hub to Venous Bloodline + Venous Catheter Hub to Arterial Catheter Hub

NO

TRANSONIC DELIVERED BLOOD FLOW RATE (QB) IS >10% LOWER THAN THE HEMODIALYSIS MACHINE'S SET BLOOD PUMP SPEED OR DELIVERY FLOW READING\*

Only proceed if both catheter lumens had blood return with treatment initiation.

Using aseptic technique, reverse the catheter configuration by reversing the bloodlines to the opposite lumens of the catheter used for the initial measurement. Document configuration.

Repeat the blood flow measurement.

NO

TRANSONIC DELIVERED BLOOD FLOW RATE (QB) IS 10% LOWER THAN THE HEMODIALYSIS MACHINE'S SET BLOOD PUMP SPEED

Carefully document measurement and catheter configurations.


Proceed to recirculation measurements with both catheter configurations.

Escalate the results of the findings to the nephrologist for possible catheter evaluation or prescription adjustment to address catheter dysfunction.

# Optimizing HD Adequacy in Catheters

## Catheter Configuration with the Transonic HD Monitor


### Step 2:



# Optimizing HD Adequacy in Catheters cont.

## Catheter Configuration with the Transonic HD Monitor For Use with Fresenius 5008 or other Hemodialysis machines that have Compensated Blood Flow Rate Capabilities

### Step 1:




# Optimizing HD Adequacy in Catheters cont.

## Catheter Configuration with the Transonic HD Monitor

For Use with Fresenius 5008 or other Hemodialysis machines that have Compensated Blood Flow Rate Capabilities cont.

Step 2:



## Catheter References

Evanson JA *et al*, "Measurement of the Delivery of Dialysis in Acute Renal Failure," *Kidney Int* 1999; 55: 1501- 508.

Sands JJ *et al*, "Difference between Delivered and Prescribed Blood Flow (QB) in Hemodialysis," *ASAIO J* 1996; 42(5): M717-719.

*Link to the 2019 KDOQI Vascular Access Guidelines*

[https://www.kidney.org/professionals/guidelines/guidelines\\_commentaries/vascular-access](https://www.kidney.org/professionals/guidelines/guidelines_commentaries/vascular-access)

*Guideline Implementation Toolkit for Monitoring and Prevention of CVC Dysfunction*

[https://www.kidney.org/sites/default/files/vait-20\\_cvc\\_complications-monitoring\\_detection\\_cvc\\_dysfunction.pdf](https://www.kidney.org/sites/default/files/vait-20_cvc_complications-monitoring_detection_cvc_dysfunction.pdf)

European Renal Association-European Dialysis and Transplant Association (ERA-EDTA), *European Best Practice Guidelines on Haemodialysis: Guideline 5.2. Nephrol Dial Transplant* 2007; 22(Suppl 2): ii99-ii100.

KHA-CARI Guideline: Vascular access – central venous catheters, arteriovenous fistulae and arteriovenous grafts.

Polkinghorne KR, Chin GK, MacGinley RJ, Owen AR, Russell C, Talaulikar G, Vale E, Lopez-Vargas PA. *Nephrology* 2013; 18(11): 701-5.

Carson, RC, Macrae, J, Kiaii, M, "Blood Pump Speed, Recirculation, and Urea Clearance in Hemodialysis Patients with Dysfunctional Catheters," *ASAIO Journal*, 2003 Abstracts, Vol 49, No 2, p 200, 2003.

Mulec, H, Henriksson, E-L, Fransson, E, Dahlberg, P, "Meticulous Medical Care Provides Minimum of Complications and Excellent Survival of Tunneled Central Venous Catheters," *ASN Abstracts* 2004 [SA-PO313]

Level, C, Lasseur, C, Chauveau, P, Bonarek, H, Perrault, L, Combe, C, "Performance of twin central venous catheters: Influence of the inversion of inlet and outlet on recirculation," *Blood Purification*, Vol 20, No2, p 182-188, 2002.

Leblanc, M, Bosc, J-Y, Vaussenat, F, Maurice, F, Leray-Moragues, H, Canaud, B, "Effective Blood Flow and Recirculation Rates in Internal Jugular Vein Twin Catheters: Measurement by Ultrasound Velocity Dilution," *American Journal of Kidney Diseases*, Vol 31, No 1, p 87-92, 1998.

Leblanc, M, Bosc, J-Y, Paganini, EP, Canaud, Advances in Renal Replacement Therapy, 1997; 4(4): 377-389.

Trerotola, SO, Kraus, MA, Shah, H, Namyslowski, J, Johnson, MS, Stecker, MS, Patel, N, "Randomized comparison of split tip Vs step tip high flow hemodialysis catheters," *JASN Abstracts*, Vol 12, p 305A, A1568, 2001.

Beathard, GA, Jefferson, VA, Carter, MJ, "Clinical Evaluation of a Subcutaneous Dialysis Access Port," *JASN Abstracts*, Vol 9, p 167A, 1998.

Trerotola, SA, Kraus, M, Gassensmith, C, Ambrosius, WT, "Randomized Study of Conventional Versus High Flow Hemodialysis Catheters," *JASN Abstracts*, Vol 9, p 185A, 1998.

Kapoian, T, Syed, ST, Sherman, A, "An Assessment of Central Vein (CV) Hemodialysis Catheter Function," *JASN Abstracts*, Vol 8, p 161A, 1997.

Leblanc, M, Bosc, F, Vaussenat, F, Leray, H, Maurice, F, Gerred, LJ, Canaud, B, "Effect Blood Flow and Recirculation Rates on Internal Jugular Vein Twin Catheters: Measurement by Ultrasound Dilution," *JASN Abstracts*, Vol 8, p 165A, 1997.

Shapiro, W, Gurevich, L, "Inadvertent Reversal of Hemodialysis Lines - A Possible Cause of Decreased Hemodialysis ((Transonic Reference #) Efficiency," *JASN Abstracts*, Vol 8, p 172A, 1997.

Brugger, J, Finch, D, "Flow Rate Comparisons for Long Term, Cuffed, Central Venous Catheters and a Novel Subcutaneous Venous Access System (VAS)," *JASN Abstracts*, Vol 8, p 154A, 1997.

Burbank, J, Finch, D, "A New Subcutaneous Vascular Access System for Hemodialysis and Apheresis: Animal Studies," *JASN Abstracts*, Vol 8, p 154A, 1997.

Sands, J, Jabayc, P, Miranda, C, "Delivered Blood Flow in Cuffed Central Venous Dialysis Catheters," *ASAIO Abstracts*, Vol 43, No 2, p 69, 1997.

Mandelbaum, AP, Fusser, M, Wiesel, M, "Efficacy of a New Double-Lumen, Unipuncture Hemodialysis Access Needle with Continuous Blood Flow - First Results," *JASN Abstracts*, Vol 10, p 211A, 1999.

Hachicha, M, Huu, TC, Bellou, L, Cannard, L, Kessler, M, "Permanent Catheter Implantation Via a Persistent Left Superior Vena Cava," *JASN Abstracts*, Vol 13, p 711A, 2002.

Dipchand, CS, Jindal, KK, Keough-Ryan, TM, Thompson, KJ, Hirsch, D, "The Relationship between URR and Recirculation in Temporary Symmetric Dialysis Catheter-Flow and Recirculation in a Swine Model," *JASN Abstracts*, Vol 14, p 242A, 2003 .

Mankus RA1, Ash SR, Sutton JM, *Comparison of blood flow rates and hydraulic resistance between the Mahurkar catheter, the Tesio twin catheter, and the Ash Split Cath,* *ASAIO J*. 1998 Sep-Oct;44(5):M532-4.



[www.transonic.com](http://www.transonic.com)

Transonic Systems Inc. is a global manufacturer of innovative biomedical flow measurement equipment. Founded in 1983, Transonic sells state-of-the-art, transit-time ultrasound devices for surgical, hemodialysis, perfusion, ECMO, and medical device testing applications, and for incorporation into leading edge medical devices..

### AMERICAS

Transonic Systems Inc.  
34 Dutch Mill Rd  
Ithaca, NY 14850  
U.S.A.  
Tel: +1 607-257-5300  
Fax: +1 607-257-7256  
support@transonic.com

### EUROPE

Transonic Europe BV  
Business Park Stein 205  
6181 MB Elstloo  
The Netherlands  
Tel: +31 43-407-7200  
Fax: +31 43-407-7201  
europe@transonic.com

### ASIA/PACIFIC

Transonic Asia Inc.  
6F-3 No 5 Hangsiang Rd  
Dayuan, Taoyuan County  
33747 Taiwan, R.O.C.  
Tel: +886 3399-5806  
Fax: +886 3399-5805  
support@transoniasia.com

### JAPAN

Nipro-Transonic Japan Inc.  
7th Floor, Maruha Building  
11-1 Matsuba-cho  
Tokorozawa City, Saitama  
359-0044 Japan  
Tel: +81 04-2946-8541  
Fax: +81 04-2946-8542  
info@transonic.jp