








Transonic Inside™ Measurement Solutions for ECMO





ECMO Experts


Transonic is proud to have been part of some of the most important ECMO developments.



ECMO Measurement Technology

Transit-time Ultrasound
 +
 Ultrasound Dilution Technology
 =
 A Uniquely Powerful Combination!


Flow/ Dilution Sensor:
 The blood line is inserted into the channel of the flow/dilution sensor body. The direction of flow is indicated by arrows. Ultrasonic beam is shown emanating from the two transducers in the



Transit-time Ultrasound Flow Measurement Methodology

Using wide-beam illumination, transducers pass ultrasonic signals back and forth, alternately intersecting flowing liquid in upstream and downstream directions. The transit time of the ultrasonic beam is decreased when traveling downstream and increased when traveling upstream.

The integrated differences between the upstream and downstream transit times over the distance of the tubing/vessel provide a measure of volume flow. The sum of the transit times yields information about liquid properties that affect ultrasound velocity and allow for a suite of ultrasound dilution measurements.



Ultrasound Dilution Technology Measurement Methodology

1. A bolus of isotonic saline is introduced into the extracorporeal circuit, diluting the blood. A corresponding drop in protein concentration reduces the liquid's ultrasound velocity.
2. Paired flow/dilution sensors (arterial and venous) attached to the extracorporeal circuit sense this decrease in velocity. Indicator dilution curves can then be produced and displayed.
3. This flow-derived measurement data allows for a variety of parameters to be calculated including recirculation, clotting indication, cardiac function measures and much more.

The Krivitski Method[®]

Pioneered by Nikolai Krivitski, PhD, DSc, the Krivitski Method[®] marries transit-time ultrasound and indicator dilution principles. It can be used in any situation involving an extracorporeal blood circuit and leverages Transonic's highly accurate tubing flow sensors. Measurements are immediate, non-invasive, and can provide a wide variety of clinically valuable data points for healthcare providers. It is commonly used today in the hemodialysis and ECMO fields.



Nikolai Krivitski

ECMO Parameters

Through the use of UDT technology leveraging dual sensors and bi-directional flows, it becomes possible to acquire many parameters beyond just volume flow:

Standard Measurements Taken During ECMO

- ✓ Mean Arterial BP (MAP)
- ✓ Oxygen Saturation
- ✓ CO₂, End-tidal CO₂
- ✓ Central Venous Pressure (CVP)
- ✓ SVO₂ (venous sat)
- ✓ Pressure pre- and post-oxygenator
- ✓ Estimated flow based on pump data (RPM's)



Additional Measurements Possible Using UDT



Recirculation



Oxygenator Blood Volume (clotting indicator)



Cardiac Output and other cardiac function parameters*



Lung function assessment*



Additional custom parameters*

* Transonic R&D work ongoing

Recirculation

During venovenous ECMO, some reinfused oxygenated blood can be withdrawn through the drainage cannula, thus not contributing to systemic oxygen delivery. Knowing the amount of recirculation helps clinicians make informed judgements relative to treatment efficacy and adjustments that can be made (e.g. cannula positioning, pump settings, tubing issues, etc.)

Oxygenator Blood Volume/Clotting Indicator

Oxygenator clotting is a major concern in medical devices circulating blood for any length of time. A quantitative indicator of clotting provides the ability to identify issues early and avoid adverse or catastrophic outcomes. Knowing OXBV brings objective data to the bedside, rather than having to rely only on subjective visual cues for decision making.

Cardiac Output and Other Cardiac Function Parameters

When managing hemodynamically unstable patients, it can be important for clinicians to recognize changes in cardiac function. Transonic's R&D team has discovered simple, non-invasive methodologies to accurately measure Cardiac Output (CO) and other cardiac function parameters. This information can contribute to recognizing potential warning signs and improving overall therapy effectiveness.

Additional Parameters

Research is ongoing with regard to many additional parameters and Transonic is committed to exploring and developing clinically meaningful measurements that can be utilized by our OEM partners.

Let's Get Started

Transonic can help you develop meaningful ECMO measurements that will set you apart from the competition. The top flow minds in the business are at your disposal – scientists, engineers, and clinical specialists. Put us to the test with any project – big or small!

Mason Caplin

Transonic Inside™

Applications Engineer



mason.caplin@transonic.com

John Haberstock

Transonic Inside™ Market

Development Manager

john.haberstock@transonic.com



References

Krivitski N, Galyanov G, Cooper D, Said MM, Rivera O, Mikesell GT, Rais-Bahrami K, "In vitro and in vivo assessment of oxygenator blood volume for the prediction of clot formation in an ECMO circuit (theory and validation).," *Perfusion* 2018; 33(15): 51-56. (Transonic Reference # ELS11317)

Ruß M, Weber-Carstens S, Skrypnikov V, Taher M, Francis RCE, Boemke W, Pickerodt PA, "Optimization of V-V ECMO circuit determined by blood recirculation measurements improved systemic oxygenation in a 10-year-old patient," 5th Annual Symposium on ARDS, Berlin, Germany, June 2019 (Transonic Reference: ELS10569A)

Said MM, Mikesell GT, Rivera O, Khodayar Rais-Bahrami K, "Precision and Accuracy of the New Transonic ELSA Monitor to Quantify Oxygenator Blood Volume (in-vivo and in-vitro studies)," 2015 PAS Annual Meeting and Eastern SPR Annual Meeting. (Transonic Reference # ELS10230V)

Said MM, Mikesell GT, Rivera O, Khodayar Rais-Bahrami K, "Influence of central hemodynamics and dual-lumen catheter positioning on recirculation in neonatal veno-venous ECMO," 2015 PAS Annual Meeting and Eastern SPR Annual Meeting. (Transonic Reference # ELS10231V)

Palmér O, Palmér K, Hultman J, Broman M, "Cannula Design and Recirculation During Venovenous Extracorporeal Membrane Oxygenation," *ASAIO J.* 2016; 62(6): 737-742. (Transonic Reference # EC11034AH)

Körver EP, Ganushchak YM, Simons AP, Donker DW, Maessen JG, Weerwind PW, "Quantification of recirculation as an adjuvant to transthoracic

echocardiography for optimization of dual-lumen extracorporeal life support," *Intensive Care Med.* 2012 38(5): 906-909. (Transonic Reference # ELS9679AH)

Clements D, Primmer J, Ryman P, Marr B, Searles B, Darling E. "Measurements of recirculation during neonatal veno-venous extracorporeal membrane oxygenation: clinical application of the ultrasound dilution technique." *J Extra Corpor Technol.* 2008;40(3):184-7. (Transonic Reference: ELS9680AH)

Walker JL, Gelfond J, Zarzabal LA, Darling E, "Calculating mixed venous saturation during veno-venous extracorporeal membrane oxygenation," *Perfusion.* 2009; 24(5): 333-9. (Transonic Reference # 7904A)

Melchior R, Darling E, Terry B, Gunst G, Searles B, "A novel method of measuring cardiac output in infants following extracorporeal procedures: preliminary validation in a swine model," *Perfusion* 2005; 20(6): 323-7. (Transonic Reference # 10234A)

Darling EM, Crowell T, Searles BE, "Use of dilutional ultrasound monitoring to detect changes in recirculation during venovenous extracorporeal membrane oxygenation in swine," *ASAIO J* 2006; 52(5): 522-4. (Transonic Reference # 7309A)

van Heijst AF, van der Staak FH, de Haan AF, Liem KD, Festen C, Geven WB, van de Bor M, "Recirculation in double lumen catheter veno-venous extracorporeal membrane oxygenation measured by an ultrasound dilution technique," *ASAIO J.* 2001; 47(4): 372-6. (Transonic Reference # HD49V)

transonic
THE MEASURE OF BETTER RESULTS.

transonic.com/oem