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
Cleaning & Sterilization of Transonic® Flowprobes

Device	All Transonic® Reusable Clinical Perivascular Flowprobes		
	PROBE TYPE	STERILIZATION METHOD	POUCH SIZE/ WRAP
	AU	ETO/Sterrad/Steris	Minimum 4x12 inch
	FMC	Steam/ETO/Sterrad/Steris	Minimum 4x12 inch
	FME	Steam/ETO/Sterrad/Steris	Minimum 4x12 inch
	FMV	Steam/ETO/Sterrad/Steris	Minimum 4x12 inch
	MR/MR-S	ETO/Sterrad/Steris	Minimum 4x12 inch
	MU	ETO/Sterrad/Steris	Minimum 4x12 inch
Warnings	Flowprobes are delicate precision instruments and should be handled carefully at all times. It is critical that the Probe connector be completely dry before use. Air dry or carefully wipe with disposable cloth/paper.		
Limitations and restrictions on reprocessing	Probes with limited reuse are programmed with the number of times they can be reused. Repeated processing has minimal effect on instruments. End of life is normally determined by wear and damage due to use. NOTE: Care must be taken when scrubbing the softer materials near the probe head to prevent damage to the silicone.		
Preparations at the point of use	Remove excess debris with disposable cloth/paper. Wipe and/or rinse with water to remove excess bio-materials.		
Containment and transportation	No particular requirements. NOTE: It is recommended that instruments are processed as soon as is reasonably practical following use. Dried-on materials are more difficult to remove.		
Cleaning preparations	No particular requirements.		
Cleaning solutions	Alkaline, Neutral or Enzymatic. Use only those cleaning agents approved by your governing regulatory agency. Use all cleaning agents according to the manufacturer's directions.		
Manual cleaning	<p>[1] Rinse excess soil from instrument (temp <30°C, 86°F).</p> <p>[2] Using detergent (e.g. Steris Prolystica 2X concentrate Neutral Detergent) and soft brush remove any visible foreign material on all probe and handle surfaces for 3 to 5 minutes. Soaking or immersion in detergent during brushing is allowed. NOTE: Excess or aggressive scrubbing of the probe neck can damage the silicone, particularly where it is sealed to the handle. NOTE: Connector surfaces may be wiped clean with solutions, but take care not to damage connector pins. If solution gets on pins, carefully wipe them dry as soon as possible.</p> <p>[3] Rinse with tap water.</p> <p>[4] Visually inspect for cleanliness and repeat cleaning if necessary.</p>		
Automatic cleaning	<p>Only use cleaning solutions which have been approved for use with an automatic washer (e.g. Steris Prolystica 2X Concentrate Alkaline Detergent).</p> <p>[1] Detergent wash minimum of 2 minutes in hot tap water</p> <p>[2] Rinse minimum of 2 minutes at 70°C</p> <p>[3] Dry minimum of 15 minutes at 80 °C</p> <p>[4] Visually inspect for cleanliness and repeat cleaning if necessary.</p> <p>NOTE: Do not exceed 90°C unless the device has the autoclave label on the connector.</p>		
Disinfection	After cleaning, all Probes must be sterilized. Additional disinfection is not required and may ultimately damage the Probe. Use only those disinfecting solutions approved by your governing regulatory agency. When disinfection is performed, follow the manufacturer's instructions applicable to the disinfection solution.		
Packaging for sterilization	<p>A polyethylene/tyvek pouch sized per the table above may be used provided it is approved by the appropriate regulatory agency for use with the desired sterilization method. Ensure that the pack is large enough to contain the instrument without stressing the seals. US customers should use a pouch that is validated and FDA cleared for the specified sterilization cycle.</p> <p>Use FDA approved sterilization wrap to cover instrument tray for Sterrad sterilization according to manufacturer's instructions.</p>		

Cleaning & Sterilization of Transonic® Flowprobes

Sterilization	<p>STERRAD</p> <p>STERRAD 100 STERRAD 100s: Short cycle STERRAD 100NX: Standard cycle STERRAD NX: Standard cycle STERRAD 200: Short cycle</p> <p>Follow the instructions for use provided with the STERRAD machine for proper sterilization processing.</p>	<p>Ethylene Oxide (ETO)</p> <p>PRECONDITIONING Humidity: 55-75% RH Temp: 38-50°C (100-122°F) Time: 12 hours</p> <p>EXPOSURE (600±50 mg/L, 3 hours) Vacuum: 0.8 ± 0.5 "HgA Sterilant gas: 100% EO Humidity: 2.4 ± 0.5 "HgA Temp: 49-54°C (120-130°F) Time: 3-3.5 hours</p>	<p>POST EXPOSURE Vacuum: 1.0 ± 0.5 "HgA</p> <p>AERATION Temp: 43-55°C (109-131°F) Time: 12 hours</p>
	<p>STERIS</p> <p>V-PRO MAX Non-lumen, cycle time = 28 min Lumen, cycle time = 60 min</p> <p>V-PRO 1 PLUS Non-lumen, cycle time = 28 min Lumen, cycle time = 60 min</p> <p>V-PRO 60 Non-lumen, cycle time = 28 min Lumen, cycle time = 60 min</p>	<p>STEAM (AUTOCLAVE)</p> <p>Only reusable Flowprobes with this symbol on the connector can be autoclaved.</p> <p>GRAVITY DISPLACEMENT STEAM STERILIZATION: 132°C for 15 minutes with 30 minutes dry time 135°C for 10 minutes with 30 minutes dry time</p> <p>DYNAMIC AIR REMOVAL STEAM STERILIZATION: 132°C for 4 minutes with 20 minutes dry time 134°C for 3 minutes with 20 minutes dry time</p>	
Inspection, maintenance and testing	<p>Inspect each Perivascular Probe for:</p> <ul style="list-style-type: none"> • A bent reflector (the reflector should be at a right angle to the Probe body). • Cracks or chips in the plastic Probe body. • Nicks in the Probe cable (if nicks are observed, do not reuse). • Damage to the silicone seal (if integrity of the silicone is compromised, do not reuse). <p> Consult the Flowmeter's Operator's Manual for testing instructions.</p>		
Storage	<p>The probe is ready for use after sterilization is complete. There are no additional storage requirements.</p>		

The instructions provided above have been validated by the manufacturer for preparing a device for re-use. It remains the responsibility of the reprocessor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility to achieve the desired result. Your reprocessing procedure should comply with local regulations.

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