

T400-Series Surgical Protocol

Sheep Carotid Artery: Chronic Blood Flow Measurement

APPLICATION BASICS

Site:	Carotid artery
Species:	Sheep
Weight:	40 kg
Duration:	Chronic
Vessel Diameter:	6 - 8 mm

PROBE

Size:	8 mm (side exit)
Reflector:	L with sliding cover
Connector:	10-pin
Cable Length:	60 cm
Catalog #:	MC-8PSS-LS-WC60-CRA10-GC

FLOWMETER

TS420 Perivascular Module

Flow Ranges Observed




Fig. 1: Instantaneous flow in the carotid artery ranged from 150 to 550 ml/min.

Application

This protocol was developed for the early validations of transit time technology. In some experiments, the carotid artery was completely exteriorized in a skin loop and a cannula was installed. The carotid artery is also a convenient location for the installation of vascular access ports and pressure catheters.

Surgical Approach




Withhold feed for 24 hours prior to surgery. Premedicate with 0.4 g glycopyrrolate IM, induce with 1 g ketamine IM and maintain anesthesia on 1.5% - 2% halothane.

Have a non sterile assistant compress the vein at the thoracic inlet so that the jugular vein may be visualized. Make a 5 cm skin incision ventral to the jugular vein as shown in Figure 2. Palpate for the pulsating artery and separate the sternomastoideus and cleidomastoideus muscle to expose the carotid sheath. Carefully incise the sheath and separate the vagus nerve from the common carotid artery. Place the L bracket of the Probe around the artery. Close the slide and secure the screw.

Appose the separated muscle layers with 2-0 non-absorbable suture as shown in the figure on the next page. Secure the Probe with sutures through the holes in the bracket and around the cable as shown in Figure 4. Create a subcutaneous tunnel to a second skin incision 5 cm from the dorsal midline of the neck. Pass the cable through the subcutaneous tunnel and exteriorize with a skin button.

Close the skin with simple interrupted sutures. Also close the stab incision and the cable exit site with simple interrupted sutures. Suture the cable to the skin for strain relief and close the incision at the cable exit site, spray the site with a topical antibiotic and place a Stent bandage over the incision. The connector is soldered on the cable during recovery.

Sheep Carotid Artery: Chronic Blood Flow Measurement Cont.



ACKNOWLEDGEMENT

Dr. Alan Dobson, Department of Physiology, NYS College of Veterinary Medicine, Cornell University, Ithaca, NY 14853

REFERENCES

- Burton, R.G. and Gorewit, R.C., "Ultrasonic Flowmeter - Uses Wide-Beam Transit Time Technique," *Medical Electronics*, 1984; 15(2): 68-73.
- Drost CJ, Dobson A, Sellers AF, Barnes RJ, Comline RS, "An Implantable Transit-Time Ultrasonic Flowmeter for Long Term Measurement of Blood Volume Flow", *Federation Proceedings*, 1984; 43(3) 538.
- Dougherty RW, "Experimental Surgery in Farm Animals," Iowa State University Press, Ames 1981.



Transonic Systems Inc. is a global manufacturer of innovative biomedical measurement equipment. Founded in 1983, Transonic sells "gold standard" transit-time ultrasound flowmeters and monitors for surgical, hemodialysis, pediatric critical care, perfusion, interventional radiology and research applications. In addition, Transonic provides pressure and pressure volume systems, laser Doppler flowmeters and telemetry systems.

AMERICAS

Transonic Systems Inc.
34 Dutch Mill Rd
Ithaca, NY 14850
U.S.A.
Tel: +1 607-257-5300
Fax: +1 607-257-7256
support@transonic.com

EUROPE

Transonic Europe B.V.
Business Park Stein 205
6181 MB Elsloo
The Netherlands
Tel: +31 43-407-7200
Fax: +31 43-407-7201
europe@transonic.com

ASIA/PACIFIC

Transonic Asia Inc.
6F-3 No 5 Hangsiang Rd
Dayuan, Taoyuan County
33747 Taiwan, R.O.C.
Tel: +886 3399-5806
Fax: +886 3399-5805
support@transonicasia.com

JAPAN

Transonic Japan Inc.
KS Bldg 201, 735-4 Kita-Akitsu
Tokorozawa Saitama
359-0038 Japan
Tel: +81 04-2946-8541
Fax: +81 04-2946-8542
info@transonic.jp