

Transonic Systems, Inc.

N-Series Probes

For use with T106 / T206 and TS410 flowmeters

*** Please include the following information on your Purchase Order ***

Contact Name: _____ Date: _____
Institution: _____ Phone: _____
Probe Size: _____ Fax: _____
Flowmeter Serial #: _____ Email: _____

Probe Calibration:	Fluid	Temp (° C)	Avg. Flow Rate	Peak Puls. Flow
Slot 1:	_____	_____	_____	_____

Definitions for above terms listed below:

Standard Calibrations:	Fluid	Temp (° C)
<i>no additional charge</i>	Blood	37
	Glycerine (40%) / Water (60%) by volume	23 or 37
	Krebs Buffer Solution	23 or 37
	Saline / Water	23 or 37

Custom Calibration:

Please Note: There is a custom fluid / temperature fee of \$75.
Sample may be required.

Flow Rates:

Average Flow Rate: _____

Peak Pulsatile Flow: _____

Customizations: *** Please Note on your Purchase Order if applicable ***

1. This probe will be used in or near a magnet (i.e. NMR, MRI...)
2. This probe will need multiple calibrations (please specify additional calibrations)

For N-Series probes used on TS410 flow modules only:

There are four calibration slots available on each probe.
One calibration slot at standard fluid / temperature is included in probe price.
Additional calibration slots are available for \$75 / additional slot.
Custom calibrations are available for \$75 / custom calibration.

Additional Calibrations:	Fluid	Temp (° C)	Avg. Flow Rate	Peak Puls. Flow
Slot 2:	_____	_____	_____	_____
Slot 3:	_____	_____	_____	_____
Slot 4:	_____	_____	_____	_____