

7,Oct,2025

Transonic Optima® Meter - 10-Year Replacement Recommendation

Important Notice for Clinicians

Transonic's Optima Flowmeters are designed with an expected service life of 10 years from the date of purchase. Beyond this point, we strongly recommend replacement to ensure continued measurement accuracy and reliability.


Over time, electronic components naturally drift and wear:

- Measurement drift can occur as internal circuits age, gradually impacting accuracy.
- Component wear (such as capacitors and backup batteries) can reduce stability and performance.
- Overall reliability decreases as devices pass their intended service life, increasing the risk of unexpected downtime.

For these reasons, after 10 years, replacement is the best way to maintain confidence in your Optima system's performance and to align with good medical device practice.

If you need assistance identifying the age of your meter or planning for replacement, please contact your Transonic Sales Representative or Customer Service staff.

Sincerely,



Miriam Tenorio
President & CEO
Transonic Systems Inc